



Full flexibility with Xplan Design Service



Clear aligners have revolutionized the orthodontic landscape in recent years. Thanks to advanced digital technology you are now able to manufacture clear aligners in-house and offer cost-efficient treatments to patients. But you don't have to go full-blown digital in one big step. In fact, the opposite is usually the best tactic.

Xplan, the design service from Dental Axess, gives you full flexibility and allows you to outsource time-consuming design steps whilst keeping control over costs and treatment quality at all times. As for manufacturing, you can choose to do this in-house or outsource to third-party manufacturers as required.

Benefits

Substantial time savings

Outsource planning and design to a professional planning service.

Case management

Track progress of cases with the Xplan dashboard and receive notifications when actions are required.

Instant changes

Make changes to the proposed treatment plan in the web-based planning tool and submit for manufacturing.

Patient communication

Increase patient acceptance using the web-based patient viewer.

Manufacturing support

Choose to manufacture in-house or outsource to third-party manufacturers. Direct 3D printer integration and automatic nesting for leading 3D printers.



Case Management Portal



Web-based Treatment Planning




Patient Viewer Tool



3D Printer Integration with automatic nesting

Clear aligner workflow with Xplan



The typical workflow starts with a digital impression of the patient using an intraoral scanner. The treatment is planned and models are designed using sophisticated CAD/CAM software before being 3D printed in-house or at a third-party manufacturer. The clear aligners are then thermoformed based on the printed models.

Treatment planning and designing is a time-consuming task and often seen as a barrier to in-house aligner manufacturing. With Xplan, the digital design service by Dental Axess, you can now outsource this step to specialists at a very competitive rate and benefit from a cost-efficient and seamless in-house clear aligner workflow. Simply upload the .STL file in the Xplan online portal and choose from a range of service options. The treatment plan is designed by a specialist and you can make further adjustments in our web-based planning tool if required.

For more information and to book a consultation with one of our specialists visit dentalaxess.com

Pricing

Case Preparation only

.STL files are prepared and teeth segmented allowing doctor to plan the treatment in-house

CHF 59.00

Xpress Aligner

.STL files are prepared and teeth 3-3 are aligned in 0-6 stages.
(Case may be upgraded to Limited Align during approval process)

CHF 79.00

Limited Aligner

.STL files are prepared and teeth 5-5 are aligned in 7-12 stages.
(Case may be upgraded to Unlimited Aligner during approval process)

CHF 99.00

Unlimited Aligner

.STL files are prepared and all teeth are aligned in 13 plus stages allowing doctor to make full adjustments during approval process.

CHF 129.00

Bracket Removal only

.STL files are based and brackets removed

Single Arch	Dual Arch
CHF 25.00	CHF 35.00

Prices exclude any applicable taxes

About us

Dental Axess was founded in 2011 upon the belief that digital dentistry gives dental professionals new possibilities to improve patients' experiences.

We enable customers to move from traditional working methods to partly or fully digital workflows by providing independent advice and customized digital solutions combined with specialist training and around-the-clock customer support.

dentalaxess.com

