


Superior scanning experience with Xcart

Make your scanning setup look professional and your every-day work easier with a cart that is specially designed for your intraoral scanner.

Xcart can be used with any intraoral scanner that uses a PC or laptop (Medit i500, 3Shape TRIOS etc.)



Modern Design

Xcart is an aesthetic and practical solution for your intraoral scanner setup. It is designed to ensure that only power and intraoral camera cables are visible, hiding all other cables. Its clean lines combined with quality materials make this cart stand out in a positive way in your dental office.



Open System

To grant you more flexibility Xcart can be used with a desktop or laptop computer without the need to buy any additional add-ons. This gives you more options in deciding which computer would be most suitable for your needs, and makes it easier to upgrade your computer later on.



Product Specifications



Base Size: 83 x 55 x 40 cm

Work Surface Size: 55 x 40 cm

PC Cabinet Size: 48 x 28 cm

Cart weight: 25 kg (excluding PC, monitor and UPS)

Display Support: VESA mount - up to 32"

Warranty: 1 year

PC + UPS

There is sufficient space to fit both PC and UPS in the cart. An optional shelf is included for laptop users.

Camera Position

The camera is positioned in a way that makes it easily accessible from your workspace, yet safely hidden, meaning you can't bump into it by mistake.

Full Aluminium Body

The Xcart body is constructed from 6mm thick aluminum; making it strong and reliable.

Options Available



Xcart Stand-Alone

No PC, monitor or battery included



Xcart ECO

Cart with Ultra PC



Xcart PRO

Cart with Ultra PC and touch screen monitor



Xcart FLEX

Cart with Ultra PC, touch screen monitor and battery

About us

Dental Axess was founded in 2011 upon the belief that digital dentistry gives dental professionals new possibilities to improve patients' experiences.

We enable customers to move from traditional working methods to partly or fully digital workflows by providing independent advice and customized digital solutions combined with specialist training and around-the-clock customer support.

dentalaxess.com

